



**Q6 EDGE® 2.0 X**  
shown in *Electric Blue*  
with iLevel® technology  
and Q-Logic 3 Advanced  
Drive Controls with  
Bluetooth®

# Q6 edge<sup>®</sup> 2.0 X

With features like standard 6.25 mph motors and advanced responsiveness, the Q6 Edge® 2.0 X is engineered to meet the performance needs of the most active user. The Q6 Edge 2.0 X accepts a complete range of seating and electronics options to deliver exceptional rehab capability at an excellent value, and is available with our iLevel® technology.

**QUANTUM®**

# Specifications

<b>DRIVE WHEELS</b>	Standard: 14" Solid Optional: 14" Pneumatic; 14" Knobby
<b>CASTER WHEELS</b>	Front: 6" Solid Rear: 6" Solid
<b>ANTI-TIP WHEELS</b>	N/A
<b>MAXIMUM SPEED</b> <sup>1</sup>	Up to 6.25 mph <sup>10</sup>
<b>GROUND CLEARANCE</b>	2.7" (battery tray)
<b>TURNING RADIUS</b> <sup>2</sup>	20.5"
<b>OVERALL LENGTH</b> <sup>2</sup>	35.5"
<b>BASE WIDTH</b>	24.4"
<b>SEATING SIZES</b>	
<b>Synergy®/TRU-Balance® 3 Static Seat</b>	<b>Width:</b> 10-22" <b>Depth:</b> 10-22"
<b>TRU-Balance® 3 Power Tilt</b>	<b>Width:</b> 12-22" <b>Depth:</b> 12-22"
<b>TRU-Balance® 3 Lift &amp; Tilt</b>	<b>Width:</b> 12-22" <b>Depth:</b> 14-22"
<b>SEAT-TO-FLOOR HEIGHTS</b>	
<b>Synergy® Seat</b> <sup>5</sup>	15.875-20.875"
<b>TRU-Balance® 3 Static Seat</b> <sup>5</sup>	16.5-18.5"
<b>TRU-Balance® 3 Power Tilt</b> <sup>5</sup>	17.5-19.5"
<b>TRU-Balance® 3 Lift &amp; Tilt</b> <sup>5</sup>	17.5-19.5"
<b>MANUAL RECLINE</b> <sup>6</sup>	No
<b>BATTERY SIZE</b> <sup>3</sup>	NF-22
<b>BATTERY WEIGHT</b> <sup>7</sup>	38 lbs.
<b>AVAILABLE ELECTRONICS</b> <sup>8</sup>	120A Q-Logic 3 EX 120A Q-Logic 3e
<b>BATTERY CHARGER</b>	8A Off-board
<b>MOTOR PACKAGES</b>	2-Pole Mid-Wheel 6 <sup>®</sup>
<b>WEIGHT CAPACITY</b> <sup>4</sup>	300 lbs.
<b>BASE WEIGHT</b>	138.6 lbs.
<b>MEDICARE CODE COVERAGE</b>	<b>Group 3:</b> Standard (K0848/K0849) Single Power (K0856) Multiple Power (K0861)

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- 7) Battery weight may vary +/- 3 lbs. based on manufacturer.
- 8) All Quantum Rehab® units can be ordered with specialty controls.
- 9) Call for availability
- 10) 6.25 mph speed available only with Q-Logic 3.

# Features

## STANDARD FEATURES:

- 6 mph programmable motors
- 8A off-board charger
- Q-Logic 3 EX Drive Controls
- LED fender lights
- USB charger port (SOF for K0849)
- Push handles (K0848/K0856/K0861 only)
- Depth-adjustable high back seat (K0849 only)

## OPTIONAL INCLUDED FEATURES:

- 4" wide flat-free knobby extended traction tires
- Center-mount foot platform (K0848/K0856/K0861 only)
- Swing-away footrests
- Attendant control (K0856/K0861 only)
- Q-Logic 3 iAccess module (K0856/K0861 only)

# Options

- iLevel® option allows 12" of seat elevation while at walking speed (3.5 mph)
- TRU-Balance® 3 Seating
- TRU-Comfort 2 Seating
- Trilogy vent tray<sup>9</sup>
- Compact vent tray<sup>9</sup>
- LED full lighting package
- Accu-Trac Advanced Tracking Technology

# Colors

SEAT BACKS for units without TRU-Balance® 3 are only available in *Back in Black*.



**QUANTUM**®

#1 FOR REHAB POWER

(US) 866-800-2002

| (Canada) 888-570-1113

| quantumrehab.com

